



*Life made easy with Glucorx AiDEX™*



## Features and Benefits:



No Calibration



No Finger Pricks



Affordable



Accurate



Real Time Notifications



Readings Every 5 Minutes

### GlucorX AiDEX™ Continuous Glucose Monitor (CGM):

- Readings every 5 minutes
- Smartphone App - no reader required
- Bluetooth connection
- CGM has an alert system
- No need to calibrate
- Can be worn on the abdomen or upper arm (we recommend the abdomen)
- Comfortable to wear
- GlucorX AiDEX™ has an IPX7 waterproof rating (Can be submerged up to 1 meter in water for 30 mins)



[www.glucorx.co.uk](http://www.glucorx.co.uk)



0800 007 5892

*Life made easy with GlucorX AiDEX™*


00162CL/V012203.22




## GlucoRx AiDEX™ Step-By-Step Guide:

- 

1. Download the GlucoRx AiDEX™ smartphone app

---
- 

2. Apply the GlucoRx AiDEX™ 14 day sensor to the arm or abdomen.

---
- 

3. Apply the GlucoRx AiDEX™ transmitter to the sensor. Transmitter lasts 4 years.

---
- 

4. Connect GlucoRx AiDEX™ app for data.



[www.glucorx.co.uk](http://www.glucorx.co.uk)



0800 007 5892

*Life made easy with GlucoRx AiDEX™*

00162CL/V012203.22



GlucorX AiDEX™ Sensor Pip Code: 419-6127

Drug Tariff GlucorX AiDEX™ pricing: Sensor: £29.95 per sensor  
Transmitter: Free

### CGM Advantages:

- Overall MARD\* - 9.08%
- Overall accuracy - 98.69%
- IPX7 waterproof rating (30 minutes at 1 meter depth)
- 14 day sensor lifetime
- Transmitter lasts 4 years
- DVLA compliant
- Bluetooth connection to smartphone
- Share data
- Compact size: 35mm x 21mm x 9mm

*\*overall MARD of 9.08% (when compared to venous blood reference)*



[www.glucorx.co.uk](http://www.glucorx.co.uk)



0800 007 5892

*Life made easy with GlucorX AiDEX™*

00162CL/V012203.22