

# Preferred product list for adults

## QUICK REFERENCE GUIDE

This guide is for first line choices for most catheter patients. The list does not preclude other products where there are specific clinical reasons that make a first line option inappropriate. For further guidance contact a representative:

**Continence Service:** 01823 346158

**ICB:** 01935 385276/384123

**LINC Medical Systems:** 01572 717515

**BD Urology:** 01293 606829

**Optimum Medical:** 07872 470445

**Teleflex:** 01494 532761

**UroFlush:** 07741 775726



## CONTINENCE SERVICE

Updated January 2026

Products	Company	Name	Size	PIP Code	Product Code	Price	Comments/Average Usage
<b>Catheters</b>  Standard catheters are used in the majority of patients.  Female catheters are only used in mobile ladies and must <b>never</b> be prescribed in males	<b>Teleflex</b>  Female 28cm  Standard 40cm  With syringe pre-filled with sterile water and empty 10ml syringe for catheter removal	<b>AQUAFLATE</b>  (PTFE Coated Latex)	28cm 12	285-6300	DP2101/12	£2.36	Indwelling catheter (short-term use – up to 28 days) Usual supply <b>1 per month</b>
			28cm 14	285-6318	DP2101/14		
			28cm 16	285-6326	DP2101/16		
			40cm 12	285-6821	DP3101/12		
			40cm 14	285-6839	DP3101/14		
			40cm 16	285-6847	DP3101/16		
		<b>SYMPACATH AQUAFLATE</b>  (First Choice – Hydrogel coated Latex)	28cm 12	285-6151	DH2101/12	£6.75	Indwelling catheter (long-term use 4 – 12 weeks) Usual supply <b>1 every 3 months</b>
			28cm 14	285-6169	DH2101/14		
			28cm 16	285-6177	DH2101/16		
			40cm 12	285-6029	DH3101/12		
			40cm 14	285-6060	DH3101/14		
			40cm 16	285-6078	DH3101/16		
		<b>BRILLANT AQUAFLATE</b>  (First choice – 100% Silicone)	28cm 12	285-5864	DA2101/12	£6.25	
			28cm 14	285-5872	DA2101/14		
			28cm 16	285-5880	DA2101/16		
			40cm 12	285-6227	DA3101/12		
			40cm 14	285-6235	DA3101/14		
			40cm 16	285-6243	DA3101/16		
<b>Catheters</b>  Standard catheters are used in the majority of patients.  Female catheters are only used in mobile ladies and must <b>never</b> be prescribed in males	<b>Optimum Medical - UGO Box</b> (First Choice: Silicone Catheter)  Catheter All Silicone Standard 40cm  Everything you require for removal and insertion of a urethral catheter.	Leg Bag short tube	40cm	415-0660	3450	£16.97	Indwelling catheter (Long Term use 4-12 weeks)  Usual supply: <b>1 every 3 months</b>
			12,14 or 16ch	415-0678	3451		
				415-0686	3452		
		Leg Bag long tube	40cm	415-0728	3456		
			12,14 or 16ch	415-0736	3457		
				415-0744	3458		

Products	Company	Name	Size	PIP Code	Product Code	Price	Comments/Average Usage
<b>Catheters</b> Open-ended Long-Term Use (up to 12 weeks)	<b>Optimum Medical</b> Includes prefilled sterile water syringe for inflation and an empty 10ml syringe for catheter removal.	<b>Open tip all-silicone</b> Standard length (1 per pack)	12ch (Fr) 14ch (Fr) 16ch (Fr) 18ch (Fr)	435-2431 435-2449 435-2456 435-2464	OE-UFC12 OE-UFC14 OE-UFC16 OE-UFC18	£5.80	1 per pack Usual supply <b>1 every 3 months.</b>
<b>Catheters Opti Tip</b>  Balloon closer to tip and drainage eyehole below the balloon	<b>LINC Medical</b> With syringe pre-filled with lubricant gel, glycerine solution and empty 10ml syringe for catheter removal	<b>Opti Tip All Silicone</b>	40cm 12 40cm 14 40cm 16 40cm 18 28cm 12 28cm 14 28cm 16 28cm 18	4045910 4045928 4045936 4045944 4045951 4045969 4045977 4045985	1007UB12OT 1007UB14OT 1007UB16OT 1007UB18OT 1007UBF12OT 1007UBF14OT 1007UBF16OT 1007UBF18OT	£11.36	Indwelling catheter (long-term use 4 – 12 weeks) Usual supply <b>1 every 3 months</b>
<b>Catheters Longer Length</b>	<b>LINC Medical</b> With syringe pre-filled with glycerine solution and empty 10ml syringe for catheter removal	<b>LINC All Silicone</b>	60cm 12 60cm 14 60cm 16	4367082 4367108 4367090	08501210L 08501410L 08501610L	£8.50	Indwelling catheter (long-term use 4 – 12 weeks) Usual supply <b>1 every 3 months</b>
<b>Catheter Valve</b>	<b>LINC Medical</b>	Lever Tap (5 per pack)	N/A	344-1524	CF1	£8.28	No more than <b>1pk 5 each month.</b>

Products	Company	Name	Size	PIP Code	Product Code	Price	Comments/Average Usage
<b>Incontinence Sheath</b>  (can remain in situ for up to 24hrs)	Optimum Medical	UGO Sheath Range  Silicone adhesive sheath (28 per box)	<u>Standard</u> 25mm 28mm 30mm 32mm 36mm 40mm	413-4110 413-4102 413-4128 413-4136 413-4144 413-4151	3100 3101 3102 3103 3104 3105	£49.31	1 Sheath/24 Hours
			<u>Pop On</u> 25mm 28mm 30mm 32mm 36mm 40mm	413-4169 413-4177 413-4185 413-4193 413-4201 413-4219	3106 3107 3108 3109 3110 3111		
<b>Catheter Maintenance Solutions</b>  <b>ONLY TO BE USED AFTER FOLLOWING THE CATHETER PATHWAY</b>  <b>FOR SHORT TERM USE</b>	UroFlush	Saline (Sodium Chloride 0.9%)	50ml 100ml	414-3186 414-3194	URO50S URO100S	£3.15	Can be used for mucus, tissue debris, blood clots eg after surgery
		Solution G (3.23% Citric Acid)	50ml 100ml	414-3228 414-3236	URO50G URO100G	£3.15	Prevents phosphate crystallisation and dissolves existing crystals formed by urease producing bacteria
		Solution R (6% Citric Acid)	100ml	414-3210	URO100R	£3.15	Dissolves severe encrustation – large, gritty crystal deposits on tip of catheter. Ideal to be used just prior to catheter removal

Products	Company	Name	Size	PIP Code	Product Code	Price	Comments/Average Usage
<b>Leg Bag x 4 Night Bag x 28</b> (can remain in situ for up to 7 days)	<b>Optimum Medical</b>	Ugo 4 weeks box  4 x leg bags 28 x night bags 4 x pairs leg bag straps 4 x catheter security clips	Short Tube Long Tube  500ml Lever Tap	403-8824 403-8840	3014 3015	£31.92 £31.92	4-week urinary drainage management system
<b>Leg Bag x 4 Sterile Night Bag with lever tap x 4</b>	<b>Optimum Medical</b>	Ugo 4 Weeks <b>continuous drainage</b> box  4 x leg bags 4 x sterile night bags 4 x pairs leg bag straps 4 x catheter security clips	500ml, Long Tube, Lever Tap	435-2472	3025	£31.92	4-week urinary continuous drainage management system
<b>Leg Bags</b>  (can remain in situ for up to 7 days)	<b>LINC Medical</b>  with fabric backing, a pair of straps, and a pair of gloves	LINC Flo  <b>Lever tap</b>	500ml, Direct 500ml, 10cm 500ml, 30cm 500ml, 60cm 750ml, 60cm	3712114 3712213 3783891 3811643 3834363	LM500SD-L LM500MD-L LM500LD-L LM500AD-L LM750AD-L	£24.29 £24.29 £24.29 £24.01 £24.30	Single-use  Suitable for attaching directly to a catheter for up to 7 days
<b>Night Bags</b>	<b>LINC Medical</b>  2 and 3L	FixStrip Scissors included to cut to required length (4 per pack)	2L SU Sterile 2L SU Non-sterile 3L Sterile 2L Tap, gloves 2L Lever, gloves 2L SU T-tap 3L SU T-tap	3054780 3054806 3945201 3406071 3712247 4131744 4131736	LM2LNS LM2LS LM3LS SLC13 SLC13-L LM2LST LM3LST	£2.32 £2.89 £5.72 £10.38 £10.38 £10.79 £16.19	Suitable for attaching directly to a leg bag
<b>Leg Bag Holder/Support</b>	<b>Optimum Medical</b>	Ugo Fix Sleeve Leg Bag Holder  (4 per pack)	<b>S</b> (24-39cm) <b>M</b> (36-55cm) <b>L</b> (40-70cm) <b>XL</b> (65-90cm) <b>XXL</b> (75-105cm)	394-0178 394-0160 394-0152 394-0145 409-6699	3005 3006 3007 3008 3027	£8.02   XXL £8.93	Leg bag holders/straps should be used to support leg bags & protect bladder neck from trauma. Change every 6 – 12 months

<b>Retaining Strap</b>	<b>LINC Medical</b>	FixStrip Scissors included to cut to required length (4 per pack)	Up to 100cm	FIXSTRIP	4049938	£10.86	To secure catheter tubing and prevent movement / urethral traction. Change every 6-12 months
------------------------	---------------------	---	-------------	----------	---------	--------	---

*'Please could nurses, who use the catheter trays, feedback if any products in the trays contents are not being used, as the NHS has a target to be carbon neutral'.*

**Feedback can be emailed to** - Caroline Taylor, Prescribing Technician, Directorate of Clinical, Medicine and Primary Care.

Email: [caroline.taylor26@nhs.net](mailto:caroline.taylor26@nhs.net)